


Laboratory Diagnosis of Tuberculosis by Sputum Microscopy



The handbook

Global edition

A publication of the Global Laboratory Initiative
a Working Group of the StopTB Partnership



Stop TB Partnership

Published by

SA Pathology
Frome Road
Adelaide
South Australia 5000
www.sapathology.sa.gov.au

All rights reserved. No part of this book may be reproduced or transmitted in any form or by any means without the written permission of the publisher. The authors assert their moral rights in the work.

ISBN: 978-1-74243-602-9

Production

Project Editor

Mark Fitz-Gerald

Technical Editors

Pawan Angra (CDC), Valentina Anisimova (KNCV),
Christopher Gilpin (WHO), Richard Lumb (SA Pathology),
Hiroko Matsumoto (JATA), Armand Van Deun (The Union)

Manuscript Editors

Mark Fitz-Gerald, Richard Lumb

Illustrations

Kerry Reid

Design

Sue Dyer Design

Contact

gli_secretariat@who.int

Laboratory Diagnosis of Tuberculosis by Sputum Microscopy

Richard Lumb
Armand Van Deun
Ivan Bastian
Mark Fitz-Gerald

The handbook
Global edition

Laboratory Diagnosis of Tuberculosis by Sputum Microscopy

Acknowledgements

The writing team wish to thank the World Health Organization (WHO) and the International Union Against Tuberculosis and Lung Disease (The Union) for recognizing that smear microscopy remains an important tool for the laboratory diagnosis of tuberculosis and in monitoring patients' response to treatment.

The Global Health Bureau, Office of Health, Infectious Disease and Nutrition (HIDN), US Agency for International Development, financially supports the development of *The Handbook: Laboratory diagnosis of tuberculosis by sputum microscopy* through TB CARE I under the terms of Agreement No. AID-OAA-A-10-00020. This publication was made possible by the generous support of the American people through the United States Agency for International Development (USAID). The contents are the responsibility of TB CARE I and do not necessarily reflect the views of USAID or the United States Government.

SA Pathology (formerly IMVS Pathology) is thanked for enabling the Tanzanian edition of *The Handbook* to be used as the template for the Global Laboratory Initiative (GLI) edition of *The Handbook*. GLI is a working group of the Stop TB Partnership.

Thanks go also to illustrator Kerry Reid and Sue Dyer Design for their long-standing contributions to earlier editions of *The Handbook*, and their enthusiastic and very professional inputs.

As we remarked in the first edition of *The Handbook*, it remains our sincere hope that the intended users of *The Handbook*, the technicians at the forefront of the international effort to contain and overcome TB, will find it useful in their daily laboratory work.

The handbook

Global edition



USAID
FROM THE AMERICAN PEOPLE

TB CARE I




SA
PATHOLOGY

	Page
Foreword	4
Introduction	6
Symbols and warnings	7
Personal safety	8
Sputum collection	10
Two specimens	10
Hospital patients	10
Safe collection	10
Pre-collection patient advice	11
How to collect a specimen	12
Specimen quality	12
Rejection criteria	12
Registration	13
Storage and transport	14
Smear preparation	18
What you need	18
Making a smear	19
Microscopy	
Method A Brightfield Microscopy	22
Method B Fluorescence Microscopy	50
Appendices	68
Specimen containers	68
Documentation	69
Abbreviations	72
Biosafety	73
Quality assurance	79
References	83
Patient information	84

Foreword

In 2011, there were an estimated 8.7 million new cases of tuberculosis (13% coinfecting with HIV). 1.4 million people died from the disease, including almost one million deaths among HIV-negative individuals and 430,000 among people who were HIV-positive. 5.8 million (67%) of these newly diagnosed cases were notified to national TB control programmes and reported to the World Health Organization. Among the 4.5 million new cases with pulmonary TB, 2.6 million (56%) had sputum smear-positive TB, and another 1.9 million were smear-negative*.



In many countries, sputum smear microscopy remains the primary tool for the laboratory diagnosis of tuberculosis. It requires simple laboratory facilities, and when performed correctly, has a role in rapidly identifying infectious cases. It has been shown conclusively that good-quality microscopy of two consecutive sputum specimens will identify the vast majority (95–98%) of smear-positive TB patients**. Moreover, microscopy can be decentralised to peripheral laboratories.

Despite its advantages sputum smear microscopy does fall short in test sensitivity, especially for certain patient groups such as those living with HIV/AIDS, and also in the laboratory diagnosis of childhood and extrapulmonary disease. New diagnostic tools endorsed by WHO (such as liquid culture, line probe assay, Xpert MTB/RIF) overcome many of the limitations of smear microscopy, especially for patients living with HIV/AIDS and those with a high likelihood of having drug-resistant TB.

WHO and The Union have previously published guidelines for sputum smear microscopy. In the decade since publication, many developments have occurred and a revised and updated text replacing both is timely.

The Handbook is a practical guide for the laboratory technician; it draws on the ideas outlined above and references best practice documents released by WHO and the GLI. *The Handbook* uses simple text and clear illustrations to assist laboratory staff in understanding the important issues involved in conducting sputum smear microscopy for the diagnosis of TB.

* WHO *Global tuberculosis report 2012* WHO/HTM/TB/2012.6

**WHO *Same-day diagnosis of tuberculosis by microscopy 2011*. WHO/HTM/TB/2011.7

The purpose of *The Handbook* is to teach laboratory technicians how to safely collect, process and examine sputum specimens for the laboratory diagnosis of tuberculosis (TB).

Sputum microscopy

Sputum smear microscopy is one of the most efficient tools for identifying people with infectious TB.

Smear-positive patients are up to ten times more likely to be infectious than are smear-negative patients.

The purpose of sputum microscopy is to:

- Diagnose people with infectious TB
- Monitor the progress of treatment
- Confirm that cure has been achieved

Consistent and accurate laboratory practice helps to save lives and improves public health.

Risk of infection

Where good laboratory practices are used, risk of infection to laboratory technicians is very low during smear preparation.

A higher risk of infection exists when collecting sputum specimens from patients.

Doctors and nurses working in TB wards and clinics where aerosols are generated have a much higher risk of becoming infected with TB.

Personal safety

When performed correctly sputum examination will not place laboratory technicians at increased risk of developing TB.



Failure to follow these instructions may harm your health or cause immediate damage to equipment



Failure to follow these instructions may affect test results, or cause equipment damage over time



Correct – the preferred way to do something



Do not do this



Wear gloves for this procedure



Wear a laboratory coat for this procedure



Wash your hands



This substance is toxic



This substance is corrosive




This substance is infectious



This substance is flammable

Risks and transmission

TB is an infectious disease. Transmission occurs when small aerosols containing acid-fast bacilli (AFB) become airborne and are inhaled. When a person coughs, sneezes, sings or vigorously exhales they produce aerosols that could be infectious if the person has pulmonary TB.



Properly trained technicians, when working correctly, have a very low risk of infection in a TB laboratory. However some activities such as talking to infected patients and collecting specimens can carry a greater risk.

Assume all specimens are infectious. Do not shake or stir samples, aerosols may be generated

Specimens may contain pathogens other than TB. When working in the laboratory do not:

- Put anything in your mouth (e.g. a pen, your fingers etc.)
- Eat, drink or smoke
- Pipette by mouth
- Lick labels and envelopes etc.
- Apply cosmetics or handle contact lenses
- Store food or drinks in the laboratory
- Wear open-toed footwear or bare feet
- Use mobile telephones in the laboratory

Personal protective equipment (PPE)

Laboratory staff must be supplied with PPE that is appropriate for the microscopy laboratory.

- You must wear protective clothing at all times in the laboratory
- You must wear gloves when handling specimens
- Do not take PPE out of the laboratory
- Store PPE separately from personal clothing





Gloves

Wear gloves for all procedures that may involve direct or accidental contact with sputum, blood, body fluids and other potentially infectious materials.

- After use, remove gloves and discard into the biohazard waste bin
- Wash your hands:
 - Immediately if contaminated by a sample
 - When you finish work
 - Before leaving the laboratory



Coats

A good laboratory coat protects your skin and clothing. It has long sleeves and fastens in the front. The laboratory is responsible for supplying and cleaning laboratory coats.

Masks



Surgical masks *are not designed to protect the wearer*, they are designed to stop the wearer spreading aerosols. Respirators are not required for performing sputum smear microscopy.

Aerosols

Good work practice minimises aerosol formation and contamination of work surfaces and equipment.

- *Separate* 'clean' activities (administration, microscopy) from 'dirty' activities (specimen reception, smear preparation, staining)
- *Never shake* a sputum specimen
- *Carefully open* specimen containers, the sample may have collected around the thread of the container
- Spread the sample onto the slide *gently* in a regular motion
- *Always air dry smears* before heat fixing
- Use disposable wooden applicator sticks or transfer loops for making smears
- *Always manage laboratory waste correctly*

Ventilation

Open doors and windows help reduce the risk of infection (see page 73 Biosafety).

Two specimens

Where External Quality Assessment (EQA) is well established, and staff are limited, two sputum specimens are recommended for the laboratory diagnosis of TB.

Specimen 1

- Collect the first specimen when the patient presents to the clinic
- Give the patient a labelled sputum container for the next morning's sputum collection

Specimen 2

- Patient collects early morning sputum and takes it to the clinic

Alternatively, microscopy of two consecutive sputum specimens, collected on the same day, may be performed.

Hospital patients

If the patient is in hospital, it is better to collect a sputum specimen each morning on two consecutive days.

Safe collection

Transmission of TB occurs because infectious droplets are released into the air when an infected patient coughs.

Collect specimens outside so that infectious droplets are diluted in an open, well-ventilated area

To reduce the possibility of laboratory staff becoming infected:

- Tell the patient to cover their mouth when coughing
- Collect sputum outside the laboratory, preferably outside the building and well away from other people

Do not collect sputum specimens in closed spaces like:

- Laboratories or wards
- Toilet cubicles
- Waiting rooms
- Reception rooms
- Any poorly ventilated area



Pre-collection and patient advice

Request for examination of biological specimen for TB

Treatment no: _____ Date of request: _____

Patient's name: _____ Date of Birth: _____ Sex: Male Female

Age (years): _____ Date of Birth: _____ Sex: Male Female

Home address: _____ Phone: _____

Reason for referral: Diagnosis, presumptive RR-TB MDR-TB: Yes No
OR Follow-up, if follow-up, month of treatment: _____

HIV infection?: Yes No Unknown
Previously treated for TB?: Yes No Unknown

Specimen type: Sputum Other (specify): _____

Tests requested: Microscopy Sputum MTB-RF Culture Drug susceptibility Line Probe Assay

Name and signature of requestor: _____

Microscopy results (to be completed in the laboratory)

Date sample collected	Specimen type	Laboratory serial number(s)	Visual appearance (Microscopist or self)	Result (check one)
				negative - A+ <input type="checkbox"/> DSHPHF <input type="checkbox"/> 1-12 AFB <input type="checkbox"/> 1-12 AFB (HPF) <input type="checkbox"/> 1-12 AFB (HPF) <input type="checkbox"/> ++ <input type="checkbox"/> +++ <input type="checkbox"/>


Examined by (Name and signature): _____

Date of result: _____

✓

Diagnosis
or
Follow up

- Check the Laboratory Request Form
- Fill in any missing details
- Tick Diagnosis or Follow-up



Patient advice


If dentures are present,
remove them and rinse
mouth with bottled water.

Tell the patient the best
specimen comes from
the lungs.
Saliva or nasal secretions
are unsuitable.

Patient Information
page 84



How to collect a specimen




Do not stand in front of the patient during collection

Instruct the patient to:

1. Relax, take time
2. Inhale deeply 2 to 3 times, breathe out hard each time
3. Cough deeply from the chest
4. Place the open container close to the mouth to collect the sputum
5. After collection screw the lid on tightly

Several attempts may be necessary to obtain a good quality specimen.

Specimen quality




Good quality specimen
Mucoid



Good quality specimen
Purulent



Good quality specimen
Blood stained



Poor quality specimens
are thin and watery
or composed largely
of bubbles



Keep the best sample


Rejection criteria

Repeat collection in the following cases:

- Broken or leaking specimen containers
- Specimen container details do not match the Laboratory Request Form
- The specimen has been collected into a fixative (e.g. formalin)
- Container unlabelled
- The specimen has been collected into tissue paper

Registration

Register the specimen before processing



Laboratory Register

Laboratory Request Form

1. Check patient details on the container match the Laboratory Request Form
2. Transfer patient details from the Laboratory Request Form to the Laboratory Register.
For follow-up specimens copy the Patient District Number to the appropriate column of the register
3. Write the Laboratory Number (LN) on the side of the specimen container
4. Write the LN on the Laboratory Request Form

For each patient, use the same LN and the numbers 1 and 2 to identify the:

- First specimen (1)
- Second specimen (2)

Saliva specimens must be reported on the Laboratory Request Form.

Where AFB microscopy or molecular testing for TB are not available and the patient cannot be referred, appropriate specimen storage and transport is required.

Storage

To preserve specimen quality:

- If microscopy or molecular testing only is requested refrigeration is not required
- Store specimens to be cultured in a refrigerator or keep as cool as possible
 - Do not freeze

What you need

- Permanent marker to write details on the side of the container
- Plastic bag for each specimen
- Transport box
- Master List of specimens
- Laboratory Request Forms

Approved secondary packaging (transport box) must:

- Be leak proof and strong
- Contain absorbent material, bench roll etc.
- Keep Laboratory Request Forms separate from sputum specimens
- Be kept out of sunlight

Packing checklist

Is the sputum container clearly labelled with:

- Patient name
- Date of collection
- Specimen number (1 or 2)

Always label the container never the lid (see page 11).

- Are Laboratory Request Forms completed correctly?
- Are Laboratory Request Forms packed separately from specimens?

Prepare a Master List that contains the details for each specimen being transported.

Ensure the Master List contains the name and address of the laboratory sending the specimens.


Check that the number of specimens equals that on the Master List.

Transport

- Follow local regulations for specimen transport
- Whilst delays, even in hot weather, will not affect test results, you should send packed specimens as soon as possible

Packing specimens

Put several layers of absorbent paper in the bottom of the shipping container.



1. Check each specimen container is labelled with
 - Patient name
 - Date
 - Specimen number (1 or 2)
2. Cross check specimens against Laboratory Request Form
 - Wrap each container in a separate plastic bag
3. Seal each container in a separate plastic bag
4. Put sealed bags into the shipping container
5. Put Laboratory Request Forms and Master List into a separate sealed bag
6. Pack shipping container to prevent movement
7. Add sealed bag containing the forms
8. Seal and address the shipping container



Keep cool
Store upright
Deliver urgently

- For Satellite health centres preparing sputum smears:
 - Sputum smears must be prepared as soon as possible after collection
 - Smears are easier and safer to transport than specimens
 - Couriers bring sputum smears to the Microscopy Centre for examination, and return the results
 - Avoid once-weekly courier collections because they will result in unacceptable delays
- Ensure each smear is clearly labelled and has a completed Specimen Request Form
- Keep in a slide box away from light, heat, dust, humidity, and insects
- The courier will bring the slide box and the Specimen Request Forms
- The courier should bring back an empty slide box from the Microscopy Centre
- Seal the slide box so that smears cannot fall out or break during transit

OR

- If a slide box is unavailable wrap each slide in toilet paper
- To prevent breakage put at least five slides in each bundle. Use unused slides if required



Add the first slide



And roll twice



Add one slide at a time



Until all slides are wrapped






Use a rubber band or tape to prevent unrolling



Sputum smears must be prepared promptly after collection or receipt.

To effectively prepare smears, you will need:

A dedicated solid bench with a non-absorbent surface that can be disinfected.



Never reuse sputum smear slides

Applicator stick

Bamboo/disposable applicator sticks are best because they:

- Separate purulent material from saliva faster
- Pick up more sputum
- Are faster, safer
- Are disposable, single use

Wire loops

Some technicians prefer wire loops because they can be reused however they:

- Are more time consuming
- May collect a smaller sample volume
- Are less efficient, must be flamed and cooled between samples



Never put more than one sputum specimen on each slide




Write the LN and 1 or 2 identifier on the frosted end of each slide using a pencil

For non frosted slides use a diamond pen or stylus



Aerosols may be generated


Do not mix purulent/bloodstained portions with saliva/mucous



Sputum specimen with purulent (A) portions within saliva (B)

Select only purulent or bloodstained portions of sputum

More AFB will be found in the purulent portions of a specimen.




Ziehl-Neelsen stained smear – purulent



Ziehl-Neelsen stained smear – saliva

Older sputum specimens still give excellent results for microscopy.

3



Smear the specimen in the centre of the slide, covering 2cm by 1cm

4




Discard the applicator stick into discard container after use, do not flame, do not reuse



Use a new clean applicator stick for each specimen



5



To clean a wire loop

- Insert loop in sand trap and rotate
- Flame the loop to red-hot and allow to cool


6

Retain all specimens until results are reported




Remove gloves and wash your hands after preparing smears

7



Air dry smears on a slide rack or flat surface

8




When dry, heat fix the smears:

- ensure the smear is facing upwards
- pass 3 times through the flame of a spirit lamp

Overheating will damage the bacilli




9



Correctly prepared smear

Stained smears resulting from poor smear preparation

Too thick/too big



X

Too thin



X

Not centred and too small



X

Multiples/confusing label



X


	Page
Staining	24
What you need	24
The Ziehl-Neelsen method	26
Examination	28
Reading smears	28
Appearance of acid-fast bacilli	30
Reporting	31
How to report	31
Summary	32
False-negatives	
– Consequences	
– Prevention	
False-positives	
– Consequences	
– Prevention	
Ziehl-Neelsen reagent preparation	33
The microscope	35
Trouble-shooting	42
– Staining	42
– Microscopy	44

Brightfield sputum smear microscopy requires simple laboratory facilities and is a much cheaper alternative to the complex and costly process of TB culture. However, to be effective staff must be trained, follow correct standard operating procedures, be provided with good quality equipment, consumables and reagents, and be part of a Quality Assured network of laboratories.

The Ziehl-Neelsen (ZN) technique has been the primary diagnostic technique for over 100 years. It is easier to learn to recognise ZN stained AFB compared with fluorescence microscopy. The detection of one AFB in a smear is sufficient to declare a positive result.




To stain smears using the Ziehl-Neelsen method you will need:



Staining rack to support slides over sink or bucket



Forceps



Burning stick



Slide rack for drying stained slides



Water



Timer



Staining bottle

Stain



1% carbol fuchsin

Decolouriser



25% H₂SO₄

Counterstain



0.1% methylene blue

You will require 2 – 3 volumes of decolouriser for each volume of stain

How to fold a filter

1



2




3



How to make a burning stick


1

Cotton wool




Piece of wire

2



Roll onto wire


3



Dip into alcohol


Light

1



Place the slides smear upwards, in LN order, on a staining rack over the sink or bucket, about a finger-width apart

Ensure the slides are level






Filter during use




3




- Heat each slide from below until steam rises, always keep the flame moving
- Stop heating when steam rises

Do not boil




2




Begin at the edges, cover each slide completely with carbol fuchsin

4




Leave the heated stain on the slides – minimum 10 minutes

A longer time will improve staining, provided the stain does not dry on the slide




5





- Gently rinse each slide with water
- Tilt each slide to drain off excess water


6



Do not splash adjacent slides





7




Add decolorising solution to the slide and leave for 3 minutes

8



- Gently rinse each slide with water
- Do not splash adjacent slides
- Tilt each slide to drain off excess water

9



Cover each slide with methylene blue for 60 seconds only

10



- Gently rinse each slide with water
- Do not splash adjacent slides
- Tilt each slide to drain off excess water

11



- Air dry away from direct sunlight
- Do not dry slides with blotting paper
- Clean back of slides with moist paper

12

Do not examine slides until they have dried



A correctly stained smear

Smears must be consistently and systematically examined to ensure a representative area of the smear is reported.

1




- Check the smear is facing upwards
- Apply one drop of immersion oil
- The drop must fall freely onto the smear so that the oil applicator does not become contaminated with TB organisms

Never allow the oil applicator to touch the slide




2

Use the 10X objective to focus the first smear, avoiding the oil drop. Scan the smear, looking for purulent or mucoid material. Where the smear is too thick, too thin, or contains epithelial cells only, move up or down to find purulent or mucoid material; continue scanning.



Inflammatory cells (high power) – look for areas like this



Avoid areas containing epithelial cells (low power)

3

Carefully rotate the 100X oil objective lens over the slide


4

Carefully adjust the fine focus until cells are sharp

Never allow the lens to touch the glass slide



6



Place the slides smear down on a clean piece of paper, leave overnight

Avoid contamination, always use a clean piece of toilet paper



8

Wipe the microscope lens gently with tissue paper to remove immersion oil after each positive slide and when you have finished examining a batch of slides (for cleaning agents see page 38)

5



Direction of traversing the stained slide

Examine at least 100 high power fields (one length) before recording a negative result
You should take approximately 5 minutes to read a negative smear




7




- Store the slides in LN order in a closed box. They will be needed for EQA
- Do not write the result on the slide
- Do not treat slides with xylene

- Viewed with an oil immersion lens, AFB are red, slender rods, sometimes with one or more granules
- Tubercle bacilli may occur singly, as V-shaped forms, or as clumps of bacilli
- Report fragments of bacilli – often seen during treatment


Typical morphological characteristics of *Mycobacterium tuberculosis*




Single bacilli



V-shaped forms



Clumps of bacilli



Bacilli fragments



Where possible, all positive smears should be reviewed by another technician

The number of AFB indicates how infectious the patient is. It is important to record exactly what you see.

1



Read the smear

2

What you see	What to report
No AFB in 100 fields	No AFB observed
1 – 9 AFB in 100 fields	Record exact number of bacilli
10 – 99 AFB in 100 fields	1+
1 – 10 AFB per field, check 50 fields	2+
More than 10 AFB per field, check 20 fields	3+

Laboratory Register

TB Laboratory Register for smear and Xpert MTB/RIF															
LN	Specimen ID	Date of birth	Patient address	Treatment unit	Smear Result (LN)	Xpert MTB/RIF Result (LN)	Resistance detected (LN)	Resistance detected (LN)	Resistance detected (LN)	Resistance detected (LN)	Resistance detected (LN)	Resistance detected (LN)	Resistance detected (LN)	Resistance detected (LN)	Remarks
757/1															

1. For diagnostic testing employing serial spots or other specimens the LN is the date of receipt of the first set of specimens
 2. *To use the LN to link a specimen
 3. * = previously treated, ** = not previously treated, *** = unknown
 4. Patients on TB treatment: indicate results of treatment at each follow-up examination is performed
 5. Xpert MTB/RIF test result reported as follows: **S** = MTB detected, if resistance not detected; **DR** = MTB detected, if resistance detected; **T** = MTB detected, if resistance indeterminate; **I** = AFB not detected; **F** = Invalid/ no result error
 6. Smear results reported as follows:
 Slight AFB
 Scanty (and report number of AFB) = 1-9 AFB per 100FP;
 ++ = 10-99 AFB per 100FP;
 +++ = 1-10 AFB per FP;
 **** = 10 AFB per FP
 7. If both MTB/RIF indeterminate result, indicate error code or 'Invalid'

4 Transfer the result to the Laboratory Register

757/1 biological specimen for TB

Date of request: _____

Age (years): _____ Date of Birth: _____ Sex: Male Female

Patient address: _____ Telephone: _____

Reason for examination: diagnosis, if diagnosis, presumptive TB, TB/MDR-TB; Yes No
 OR follow-up, if follow-up, month of treatment: _____

HIV infection? Yes No Unknown
 Previously treated for TB? Yes No Unknown

Specimen type: Sputum Other (specify): _____
 Test(s) requested: Microscopy Xpert MTB/RIF Culture Drug susceptibility Line Probe Assay

Name and signature of requester: _____

Microscopy results (to be completed in the laboratory)

Date sample collected (day/month/year)	Specimen type	Lab. serial no.	Result (check one)
757/1			negative: <input type="checkbox"/> 0 AFB <input type="checkbox"/> 1-9 AFB <input type="checkbox"/> 10-99 AFB <input type="checkbox"/> 1-10 AFB/FP <input type="checkbox"/> 10 AFB/FP

Examined by (Name and signature): _____ Date of result: _____

3 Record results

5 Date/Sign → 6 Return to doctor or clinic

Laboratory Request Form

1. Use the LN to find the correct patient Request Form
2. Read the smear
3. Immediately record the result on the Request Form
4. Transfer the result to the Laboratory Register
use red pen for positive results
5. Date and sign the Laboratory Request Form
6. Return the completed Laboratory Request Form to the Doctor or Clinic



Do not give results to the patients as lost reports may delay treatment
Do not write the results on the slide as they are needed for EQA checking



False-negative means reported negative but truly smear-positive

Be accurate and consistent in all your work, lives depend on you

False-negatives – consequences

- Patients with TB may not be treated resulting in on-going disease, disease transmission, or death

Prevention

- Label sputum containers, slides and laboratory forms accurately
- The specimen must contain sputum not saliva
- Select purulent material to make the smear
- Smear preparation – centred, spread evenly, 2cm x 1cm in size
- Use good quality basic fuchsin powder and reagents
- Heat carbol fuchsin until steaming
- Do not boil during fixation
- Stain with carbol fuchsin – minimum 10 minutes
- Do not overheat the carbol fuchsin
- Decolourise until no more carbol fuchsin is released, maximum 3 minutes
- Counterstain – maximum 60 seconds
- Keep the microscope well maintained and the lenses clean
- Perform regular QC on stains and reagents
- Check the slide LN matches the Laboratory Request Form before recording the result



False-positive means reported positive but truly smear-negative

Don't rush – examine at least 100 fields in one length before reporting 'No AFB observed'

False-positives – consequences

- Patients are treated or retreated unnecessarily
- Medications will be wasted

Prevention

- Ensure laboratory staff can reliably recognise acid-fast bacilli
- Label sputum containers, slides and laboratory forms accurately
- Always use new unscratched slides
- Use bamboo/wooden sticks once only
- Do not allow carbol fuchsin to dry on the smear
- Decolourise adequately
- The oil applicator must not touch the slide
- Keep the microscope well maintained, the lenses clean, store appropriately
- Perform regular QC of stains and reagents
- Check the slide LN matches the Laboratory Request Form before recording the result

Stain

Carbol fuchsin – 1.0%

Grade


Basic fuchsin powder	10g	Certified
Ethanol (or methanol)	100ml	Technical
Phenol crystals*	50g* <i>use colourless not tinted crystals</i>	Analytical
Distilled water	900ml	

**Phenol crystals and vapour are corrosive, toxic and may cause burns
Use care, prepare in a well ventilated area**

Preparation

1. Add 100ml of ethanol (or methanol) to a one litre glass flask
2. Add 50g of phenol crystals and dissolve
3. Add 10g of basic fuchsin powder
4. Mix well until dissolved
5. Add distilled water to make one litre
6. Label the bottle – **“1% carbol fuchsin”, date and initial**
8. Store in a dark bottle in a cupboard at room temperature (expiry 12 months)

Wash your hands after preparing reagents



Perform a Quality Control check and record results in the QA log book
Filter solution at time of use



Decolourising solution

Always add the acid to ethanol or water. Solutions will generate heat.

25% H₂SO₄

Grade

Concentrated sulphuric acid (H ₂ SO ₄)	250ml	Technical
Distilled water	750ml	

Preparation

- Carefully add the H₂SO₄ to the water
- Label the bottle **"25% H₂SO₄", date and initial**
- Store in a dark bottle in a cupboard at room temperature (expiry 12 months)

3% HCl in ethanol (acid alcohol)

Grade

Fuming hydrochloric acid (HCl)	30ml	Technical
95% ethanol	970ml	Technical

Preparation

- Carefully add the HCl to the ethanol
- Label the bottle **"3% HCl in ethanol", date and initial**
- Store in a dark bottle in a cupboard at room temperature (expiry 12 months)

6% HCl

Grade

Fuming hydrochloric acid (HCl)	60ml	Technical
Distilled water	940ml	

Preparation

- Carefully add the HCl to the water
- Label the bottle **"6% HCl", date and initial**
- Store in a dark bottle in a cupboard at room temperature (expiry 12 months)

Counterstain

0.1% methylene blue

Methylene blue chloride	1.0g
Distilled water	1000ml

Preparation


- Dissolve the methylene blue chloride in distilled water
- Label the bottle – **"0.1% methylene blue", date and initial**
- Store in a dark bottle in a cupboard at room temperature (expiry 12 months)

Perform a Quality Control check and record results in the QA log book.



The microscopy area should be:

- Free from dust
- On a stable level platform
- Away from centrifuges and refrigerators
- Away from water, sinks or chemicals to avoid splashes or spills
- Ergonomically correct work position (see page 77)



Setting up the microscope



For binocular microscopes with pre-centred and fixed condensers:

1. Rotate the nose-piece to the 10X objective
2. Set the variable voltage regulator to minimum
3. Turn the power on
4. Slowly adjust until the desired light intensity is reached
5. Place a stained slide onto the stage
6. Bring the smear into focus with the coarse and fine-adjustment knobs



Always use the focusing adjustment knobs to lower the stage away from the lens

7. Adjust the interpupillary distance until the right and left images merge



8. Focus the image with the right eye by looking into the right eye-piece and adjusting with the fine focus knob
9. Focus the image with the left eye by looking into the left eye piece and turning the dioptre ring
10. Open the condenser iris diaphragm so that the field is evenly lit
11. Place one drop of immersion oil onto the smear and rotate the 100X objective into it

12. Focus using the fine adjustment knob
13. Use the variable voltage regulator to achieve a comfortable illumination
14. Once the smear has been read, rotate the 100X objective away, locate the 10X objective over the slide, and then remove the slide
15. When finished, reset the voltage regulator to a minimum, and turn the power off
16. At the end of each day, use fine tissue paper to carefully remove immersion oil from the 100X lens, do not use gauze. Cover the microscope, or put it in the microscope box or return to the humidity controlled cupboard

Do's and Don'ts

- The 100X objective is the only lens requiring immersion oil
- Keep immersion oil away from other lenses
- Immersion oil must have medium viscosity and a refractive index (RI) greater than 1.5. Any synthetic, non-drying oil with an RI > 1.5 is suitable (refer to manufacturer's instructions)
- Do not use cedar wood oil as it leaves a sticky residue on the lens



Never use cedar wood oil diluted with xylene instead of immersion oil, as it will quickly destroy the lens

Immersion oil – a simple test

Good immersion oil



Poor immersion oil



**Maintenance****Do not use xylene****Cleaning lenses****Some cleaning agents will damage lenses over time – for daily cleaning use tissue paper**

Cleaning Agent	Long term use	Infrequent use
Manufacturer's recommendation	✓	✓
Ethyl ether/alcohol (80/20)	✓	✓
Alcohol	✗	✓
Benzene/petrol	✗	✓
Acetone/ketones	✗	✓
Xylene	✗	✗

- Never use xylene to clean any part of a microscope
- Remove dust and sand from dry lenses before using cleaning fluid
- When ever possible use the cleaning fluid recommended by the manufacturer
- Use a minimum amount of cleaning fluid, never dip a lens into cleaning fluid
- Fine tissue paper is best for cleaning optical surfaces as it does not scratch the lens
- Alternatively use fine quality toilet paper
- Do not use ordinary paper, or cotton wool or gauze to clean lenses
- Keep the microscope covered when not in use
- Keep the eye-pieces in place
- Fungus or dust may enter through holes where objectives in the nose-piece are missing



Cover holes from missing objectives


- If the image appears hazy with black dots, check for dust or dirt on the lenses (eye-pieces, objectives, condenser and illuminator lens). If:
 - The black dot moves when the eye-piece is rotated, then the dust is on the eye-piece
 - The black dot moves when the slide is moved, then it is on the slide
 - These two are ruled out, then assume the dust is on the objective (if inside the objective, it appears as dots; if on the outside, then as a hazy image)
- Dust can be removed using a camel-hair/artist brush or by blowing over the lens with an air brush



A simple air brush made using a Pasteur pipette and rubber bulb

Light source


- Never touch the glass bulb surface as skin oils will burn, reducing light intensity
- Use paper to hold the bulb when inserting into the microscope



Use a tissue, do not touch the bulb with your fingers


Mechanical parts

- Never disassemble the microscope – send to a specialist technician

Fungal growth

Fungus growing inside the eyepiece tube




- Fungus growth on the lenses, the eye-piece tube and prisms causes the microscope image to become hazy and unclear
- To check for fungus turn the microscope on:
 - Rotate the 10x objective into the light path
 - Take out both eyepieces, look down the eyepiece tubes for fungus
- To prevent fungal growth, the microscope should be kept in a warm cupboard or box. A cupboard with a tightly fitting door, heated by a light globe (maximum 25W), works well
 - Always leave the cupboard light on, even when the microscope is not in the cupboard



Warming box for microscope storage
Lamp 25W maximum



- If proper storage is not available, keep the microscope in the shade and with good air circulation

Correctly stained slides



Problem

Smear too pink



Cause

Insufficient decolourisation
Acid concentration very low,
or applied for too short a time

Carbol fuchsin (CF) has dried
on smear

Smear too thick

Remedy

Decolourise for longer

For commercial reagents, check
with NTP

For in-house reagents, recheck
stain preparation and QC results

Check smears are level over sink
Add sufficient CF

Prepare new smear

When correctly stained this slide
looks like **A** above

Problem

Pale acid-fast bacilli



Cause

CF prepared from poor quality reagents

CF insufficiently heated

CF staining time less than 10 minutes

Smear overheated during preparation or staining

CF reagent has expired or stored in direct sunlight

Remedy

Use reagents from reputable manufacturer

For in-house reagents, recheck preparation and QC results

Heat CF to steaming


Stain for a minimum of 10 minutes

Pass over flame 3 times, 1-2 seconds each time
Stop heating when CF steams

Replace reagent
Store stain bottle in the dark

Problem

Counterstain too dark



Cause

Excessive counterstaining time

Inadequate washing step after counterstaining

Methylene blue concentration too strong

Smear too thick

Remedy

Do not exceed 60 seconds

Extend washing step

For commercial reagents, check with NTP

For in-house reagents, recheck preparation and QC results

Prepare new smear

Problem

Deposit on slide



Cause

Stains not filtered

Soot deposit on underside of smear

Remedy

Filter stains

Clean with a moist tissue paper

Problem**Light flickers or does not turn on****Cause**

Loose plug or connection

Loose light bulb

Dirty bulb contacts

Erratic voltage supply

Faulty on-off switch

Fuse blown or transformer blown

Discoloured bulb/burnt out

Remedy

Check wall sockets, transformer, power supply

Reinstall the bulb – Do not touch bulb with fingers

Clean contacts with 70% alcohol and retry or replace bulb

Use a voltage stabiliser

Replace the switch

Replace the fuse

Replace the bulb – Do not touch bulb with fingers

Problem**Uneven illumination****Cause**

Field of view partially blocked

Iris diaphragm is almost closed or condenser is not aligned

Dirty lenses

Heavy fungal growth on lenses


Remedy

Rotate the nose-piece until it clicks into position

Recalibrate microscope

Gently wipe the lenses with lens paper/soft cloth. If the trouble persists clean with lens paper soaked in the recommended lens cleaning fluid (see page 38)

Clean the lens using lens cleaning fluid as recommended by the manufacturer

Problem**Excessive image contrast****Cause**


Iris diaphragm is almost closed

Remedy

Open diaphragm

Problem	Cause	Remedy
Unclear image with glare	Iris diaphragm too far open	Close the iris diaphragm to make the opening smaller
Specimen focused at 10x but not at higher magnification	Slide upside down	Turn it over
Specimen goes out of focus more than usual at high magnification	Slide is not flat on the stage	Clean the stage and underside of slide
Mechanical stage cannot be raised	Lock set too low	Adjust to proper height and lock
Mechanical stage is not moving, too stiff or does not stay up	Poor tension adjustment on the mechanical stage Solidified lubricants	Adjust tension with tension adjustment knob (if present) Microscope requires service

✓ **Good**



Problem

Oil immersion objective does not give a clear image



Cause

Insufficient oil on slide

Light source or condenser collector lens dirty

Poor quality immersion oil (low refractive index)

Surface of the lens is dirty

Water on slide

Bubbles in immersion oil

Oil inside lens

Remedy

Add immersion oil

Clean using lens paper and cleaning fluid

Use quality immersion oil (see page 37)


Clean lens with tissue paper
If oil/fungus inside the objective, replace lens

Air dry slides

Remove oil from slide and carefully reapply oil

Clean or replace lens

✓ Good



Problem

Dust/dirt visible in the field of view



Cause

- Dust on the collector lens of the light source
- Dust on the top-most lens of the condenser
- Dust on the eye-piece

Remedy

- Clean all surfaces
- Clean all condenser surfaces
- Clean all surfaces

Problem

Cracked objective lens



Cause

- Lens has been dropped
- Lens forced into slide or stage

Remedy

- Replace lens
- Replace lens

Problem

Regular or semi regular crescent shapes that maybe confused for AFBs

**Cause**

The glass slide is scratched

Remedy

Learn to recognise glass artefacts

Problem

Headaches/incomplete binocular vision

Cause

Improper adjustment of interpupillary distance

Dioptre adjustment was not done

Eye-pieces are not matched

Remedy

Adjust the interpupillary distance

Adjust dioptre settings

Use matched eye-pieces

Problem

Fuse blows frequently

Cause

Fuse incorrectly rated

Unstable line voltage

Remedy

Replace with correctly rated fuse

Use voltage protection device

Method | **Fluorescence
Microscopy**

B

Method **B** | Fluorescence Microscopy

	Page
Staining	52
What you need	52
Auramine method	54
Bulk staining	56
Examination	58
Reading smears	58
Appearance of acid-fast bacilli	59
Reporting	61
How to report	61
Summary	62
False-negatives	
– Consequences	
– Prevention	
False-positives	
– Consequences	
– Prevention	
Auramine reagent preparation	63
Trouble-shooting	66
– Staining	66

In 2011, WHO released a new policy on Light Emitting Diode (LED) based Fluorescent Microscopy (FM) for diagnosing TB. FM is equally accurate, at least 10% more sensitive and has qualitative, operational, cost and workload advantages for all laboratories performing sputum smear microscopy. WHO recommended a phased approach to change from brightfield microscopy to LED-based FM across the microscopy network.

LED FM offers considerable advantages over conventional FM, which requires a darkened room to read smears. Conventional FM relies on expensive mercury vapour lamps that have a limited life span, generate large amounts of heat, and are a safety hazard if broken.

For a laboratory with a high workload, bulk staining is an acceptable option and protocols are described on page 56-57.

Reporting

Due to an historical inaccuracy, the FM reporting scale for positive smears has been revised because the actual field observed is larger than previously calculated.

Low scanty positives, 1-4 AFB in one length at 200x magnification, or 1-2 in one length at 400x magnification should be confirmed by:

- viewing additional fields
- having another technician check the AFB morphology or
- collecting another sputum sample

Confirmation of FM low-positive smears by re-staining with ZN should not be done.

Quality control

AFB in FM-stained smears fade rapidly; for FM re-stain all smears. Auramine reagent must be prepared as 10X concentrated stock that keeps well for 12 months. Diluted staining solution may deteriorate within a few months, and should be prepared monthly from stock.

Introducing LED FM methods


The switch to LED FM should be carefully phased in at country level, with LED technology that meets WHO specifications. Countries using LED microscopy should retrain laboratory staff with strong emphasis on practical training of longer duration. EQA should be introduced for individual laboratories; technique validated for the network as a whole, and the effect on TB case detection rates and treatment outcomes monitored.

Staining solutions can deteriorate quickly – the solution becomes lighter





To stain smears using the Auramine method you will need:



Staining rack to support slides over sink or bucket



Forceps



Water



Slide rack for drying stained slides



Staining bottle



Timer

Stain



0.1% auramine

Decolouriser



0.5% acid-alcohol

Counterstain



0.3% methylene blue
or 0.5% potassium permanganate

You will require 1 – 2 volumes of decolouriser for each volume of stain

How to fold a filter

1







2




3






7



- Do not splash other slides
- Tilt each slide to drain off excess water

9




Gently rinse each slide with water

11



- Air dry away from direct sunlight
- Do not dry slides using blotting paper


8



Cover each slide with methylene blue for 1 minute
If permanganate is used, time is critical because a longer time may quench the acid-fast bacilli fluorescence



10



- Do not splash other slides
- Tilt each slide to drain off excess water

12

Do not examine slides until they have dried



A correctly stained smear

Bulk staining

Consider this method when workload exceeds 10 smears per day.


What you need



Timer




Slide basket



Water




Four ~600 ml glass containers able to hold the slide basket



- Stain 0.1% auramine
- Decolouriser 0.5% acid-alcohol
- Counterstain 0.3% methylene blue



Place slides in LN order, facing one direction, in a slide basket




Place the slide basket into auramine, ensure slides are covered leave for 20 minutes minimum

Do not heat



Fill the wash container with water



Remove slide basket from auramine, place into water.

Gently move the slide basket up and down (± 1 cm) 2-3 times (gentle agitation)



5



Remove from water place into decolouriser for 2 minutes – gentle agitation

7




Remove the slide basket from the decolourising solution and place into the water; gentle agitation

9

Discard the water from the container: refill and discard twice

11




Remove from water, tilt to drain

6

Rinse the wash container. Discard and refill twice

8



Remove the slide basket from the water, place into counterstain solution for 1 minute, ensure slides are covered

10

Remove the slide basket from the counterstain solution and place into the water; gentle agitation

12



Air dry slides away from direct sunlight



A correctly stained smear



Keep stained smears in the dark using a slide box or folder as fluorescence quickly fades when exposed to light

Read the smears on the same day they were stained.

AFB are stained bright yellow against a dark background, but with some filter systems they will appear green.

Use the 20X objective to scan the smear and the 40X objective for confirming suspicious objects.

Smears must be examined in a consistent way to ensure a representative area of the smear is reported. At least one length of the smear must be examined before reporting a negative result.



When the smear has been read, store the slides immediately in LN order in a closed box, as they will be needed for EQA.

Do not write the result on the slide.




Do not restrain scanty smear positives with ZN

The typical appearance of AFB is a long, slender, slightly curved rod, but variable in shape and staining intensity


They may be uniformly stained or may contain one or more gaps, or have a granular appearance. AFB occur singly, in small groups containing a few bacilli, or more rarely, as large clumps.




Long and slender



Slightly curved




Uniformly stained




One or more gaps




Granular



Single AFB



Small groups



Large clumps

Stained smears may contain fluorescing artefacts which do not have a typical bacillary shape, and sometimes also a different colour.



Non-fluorescing yellow or green coloured bacillary shapes should not be accepted as AFB

The number of AFB indicates how infectious the patient is. It is important to record exactly what you see.

What you see (200x)	What you see (400x)	What to report
No AFB in one length	No AFB in one length	No AFB observed
1-4 AFB in one length	1-2 AFB in one length	Confirmation required*
5-49 AFB in one length	3-24 AFB in one length	Scanty
3-24 AFB in one field	1-6 AFB in one field	1+
25-250 AFB in one field	7-60 AFB in one field	2+
>250 AFB in one field	>60 AFB in one field	3+

* Confirmation required by another technician or prepare another smear, stain and read

1



Read the smear

2

Laboratory Register

TB Laboratory Register for smear and Xpert MTB/RIF														
LN	Specimen	Age	Sex	DOB	Address	Treatment unit	MTB and TB Register No.	MTB infection (1,2,3)	MTB previously treated (4,5)	Resistance previously tested (6,7)	Specimen type	Xpert	Result	Remarks
757/1														

1 For diagnosis: looking for mycobacterium species or other specimens this is the date of receipt of the first set of specimens
 2 For Yes, No, No, Unknown
 3 = previously treated; No = not previously treated; Date = unknown
 4 Patient on TB treatment, indicate months of treatment as per TB follow-up examination is performed
 5 Xpert MTB/RIF test result reported as follows: **+** MTB detected, if resistance not detected; **+** MTB detected, if resistance detected; **+** MTB detected, if resistance indeterminate
 6 Smear results reported as follows: **+** No AFB; **+** Scanty (1-4 AFB per 100 FFF); **++** 5-49 AFB per 100 FFF; **+++** 50-250 AFB per 100 FFF; **++++** > 250 AFB per 100 FFF
 7 If Xpert MTB/RIF indeterminate result, indicate error code or 'invalid'

757/1 biological specimen for TB

Date of request: _____

Age (years): _____ Date of Birth: _____ Sex: Male Female

Patient address: _____ Telephone: _____

Reason for examination: Diagnosis, if diagnosis, presumptive RR-TB/MDR-TB; Yes No
 OR Follow-up, if follow-up, month of treatment: _____

HIV infection?: Yes No Unknown
 Previously treated for TB?: Yes No Unknown

Specimen type: Sputum Other (specify): _____

Tests requested: Microscopy Xpert MTB/RIF Culture Drug susceptibility Line Probe Assay

Name and signature of requester: _____

Microscopy results (to be completed in the laboratory)

Date sample collected (day/month/year)	Specimen type	Microscopy result	Result check one
757/1			

Examined by (Name and signature): _____ Date/Sign _____

Date of result: _____

4

Transfer the result to the Laboratory Register

3

Record results

5

Date/Sign

6

Return to doctor or clinic

Laboratory Request Form

1. Use the LN to find the correct patient Request Form
2. Read the smear
3. Immediately record the result on the Request Form
4. Transfer the result to the Laboratory Register
use red pen for positive results
5. Date and sign the Laboratory Request Form
6. Return the completed Laboratory Request Form to the Doctor or Clinic



Do not give results to the patients as lost reports may delay treatment
Do not write the results on the slide as they are needed for EQA checking



False-negative means reported negative but truly smear-positive

Be accurate and consistent in all your work, lives depend on you

False-negatives – consequences

- Patients with TB may not be treated resulting in on-going disease, disease transmission, or death

Prevention

- Label sputum containers, slides and laboratory forms accurately
- The specimen must contain sputum not saliva
- Select purulent material to make the smear
- Smear preparation – centred, not too thick or too small
- Use auramine solution as fresh as possible; do not prepare large quantities
- Stain with auramine – minimum 20 minutes
- Decolourise for 1-2 minutes only
- Counterstain – maximum 1 minute
- Read smears as soon as possible and keep them protected from light
- Keep the microscope well maintained and the lenses clean
- Perform QC – use positive controls every day to check staining procedure and microscope function
- Check the slide LN matches the Laboratory Request Form before recording the result



False-positive means reported positive but truly smear-negative

Don't rush – examine at least one length of a smear before recording a negative result

False-positives – consequences

- Patients are treated or retreated unnecessarily
- Medications will be wasted

Prevention

- Ensure laboratory technicians can reliably recognise acid-fast bacilli
- Label sputum containers, slides and laboratory forms accurately
- Always use new unscratched slides
- Use bamboo/wooden sticks once only
- Filter auramine staining solution during use
- Do not allow auramine to dry on the smear
- Decolourise adequately
- Keep the microscope well maintained, the lenses clean, store appropriately
- Perform QC – use positive controls every day
- Check the slide LN matches the Laboratory Request Form before recording the result



Stain

Auramine is a potential cancer causing agent – always wear gloves and clean any spills immediately

Phenol crystals and vapour are corrosive, toxic, and may cause burns; avoid contact with skin and mucous membranes, prepare in a well ventilated area

0.1% Auramine

Grade

Auramine	10.0g	Certified
Ethanol (denatured) or methanol	1000ml	Technical
Phenol crystals*	30g* <i>use colourless not tinted crystals</i>	Analytical
Distilled water	900ml	

Preparation

To ensure solutions are fresh, laboratories examining low numbers of smears should prepare smaller volumes.

Solution A

1. Add 1000ml of ethanol (or methanol) to a one-litre glass flask
 2. Add 10.0g of auramine powder, mix until dissolved completely
- Do not use heat since this can inactivate the auramine**
3. Label **"1.0% auramine in alcohol", date and initial**
 4. Store in a dark bottle in a cupboard at room temperature (expiry 12 months)



Solution B

1. Dissolve 30g of phenol crystals in 900ml distilled water, mix
2. Label the bottle **"3% phenolic solution for auramine", date and initial**
3. Store in a dark bottle in a cupboard at room temperature (expiry 12 months)

Preparation of 0.1% auramine solution

1. Add 50ml of solution A (1% auramine in alcohol) to a 500ml dark glass bottle
2. Add 450ml of solution B (phenolic solution for auramine) and mix
3. Label the bottle **"0.1% auramine", date and initial**
4. Store in a cupboard at room temperature (expiry 2 months)



Filter auramine solution when applying to smears or filling bulk staining containers.



Wash your hands after preparing reagents

Perform a Quality Control check and record results in the QA log book.



Correctly prepared auramine is a rich golden colour – discard if pale



Decolouriser

Always add the acid to ethanol. Solutions will generate heat

0.5% acid-alcohol

Grade

Fuming hydrochloric acid	5ml	Technical
Ethanol (denatured) or methanol	1000ml	Technical

Preparation

- Carefully add the hydrochloric acid to the alcohol
- Label the bottle **"0.5% acid alcohol", date and initial**
- Store in a dark bottle in a cupboard at room temperature (expiry 12 months)
- Perform QC and record in the QA log book

1% HCl in 10% alcohol in water

Grade

Fuming hydrochloric acid	10ml	Technical
Ethanol (denatured) or methanol	100ml	Technical
Distilled water	890ml	

Preparation

- Carefully add the alcohol (or methanol) to the distilled water
- Carefully add the hydrochloric acid to the 10% alcohol (or methanol) in water
- Label the bottle **"1% HCl in 10% alcohol in water"**, date and initial
- Store in a dark bottle in a cupboard at room temperature (expiry 12 months)
- Perform QC and record in the QA log book

Counterstain

Two counterstains are described; 0.3% methylene blue (preferred) is a true counterstain, whilst 0.5% potassium permanganate acts as a quenching agent.

The choice of counterstaining solution depends on the microscope system used: permanganate produces a very dark background in some systems, making it hard to keep focus. If this occurs, then 0.3% methylene blue is a better choice counterstain, although there is slightly less contrast.

0.3% methylene blue

Grade

Methylene blue	3.0g	Analytical
Distilled water	1000ml	

Preparation

- Add the methylene blue to the distilled water
- Label the bottle **"0.3% methylene blue", date and initial**
- Store in a dark bottle in a cupboard at room temperature (expiry 12 months)
- Perform QC and record in the QA log book





Potassium permanganate is a powerful oxidising agent and may cause burns

0.5% potassium permanganate

Grade

Potassium permanganate	5.0g	Technical
Distilled water	1000ml	




Preparation


1. Add the potassium permanganate to the distilled water
2. Label the bottle **“0.5% potassium permanganate”, date and initial**
3. Store in a dark bottle in a cupboard at room temperature (expiry 12 months)
4. Perform QC and record in the QA log book






The solution should be bright purple; if it is brick-red in colour it is oxidised discard it – rinse the bottle before refilling



Wash your hands after preparing reagents

Problem	Cause	Remedy
Too much fluorescence 	Insufficient decolourisation	Check decolourisation time
	Counterstain too weak or no alcohol	Prepare new reagent
	Auramine has dried on the smear	Check smears are level over sink Add sufficient stain
	Auramine not filtered	Filter auramine at time of use
	Smear too thick	Prepare new smear
	<i>Do not heat during staining</i>	


Problem	Cause	Remedy
Pale acid-fast bacilli 	Auramine has expired or stored in direct sunlight	Replace reagent Store bottle in the dark
	Auramine <0.1%	Recheck stain preparation and QC results
	Staining time <20 minutes	Stain for at least 20 minutes
	Smear overheated during fixation step	Pass smear through flame 3 times, 1-2 seconds each time
	Overdecolourised	Do not exceed the maximum time (1-2 minutes only)
	Stained smears exposed to daylight	Keep slides in the dark using slide box or similar Read smears as soon as possible
	Smear too thick	Prepare new smear

Problem	Cause	Remedy
Background too dark 	Counterstained too long (or) Decolourised too long Inadequate washing step after counterstaining Counterstain concentration too strong Smear too thick	Do not exceed 1 minute Do not exceed 2 minutes Extend washing step Ensure washing step is complete Recheck stain preparation and QC results Prepare new smear



For microscope problems refer to manufacturer's instructions

Ideal specimen container



Laboratory Request Form

This example shows the type of information required on a Specimen Request Form.

Request for examination of biological specimen for TB

Treatment Unit: _____ Date of request: _____

Patient name: _____

Age (years): _____ Date of Birth: _____ Sex: Male Female

Patient address: _____

_____ Telephone: _____

Reason for examination:

 Diagnosis. If diagnosis, presumptive RR-TB/MDR-TB?: Yes NoOR Follow-up. If follow-up, month of treatment: _____HIV infection?: Yes No UnknownPreviously treated for TB?: Yes No UnknownSpecimen type: Sputum Other (specify): _____Test(s) requested: Microscopy Xpert MTB/RIF
 Culture Drug susceptibility Line Probe Assay

Name and signature of requestor: _____

Microscopy results (to be completed in the laboratory)

Date sample collected (filled by requestor)	Specimen type	Laboratory serial number(s)	Visual appearance (blood-stained, mucopurulent or saliva)	Result (check one)				
				Negative (0 AFB / 100HPF)	1-9 / 100HPF (scanty; report number of AFB)	+	++	+++

Examined by (Name and signature): _____

Date of result: _____

Laboratory Register

This example shows the type of information required on a Laboratory Register.

TB Laboratory Register for smear and Xpert MTB/RIF

Lab. serial No.	Date specimen received ¹	Patient name	Sex M/F	Age Date of birth	Patient address	Treatment unit	BMU and TB Register No.	HIV infection (Y/ N/ Unk) ²	Patient previously treated for TB ³	Examination type (Tick one option)	Examination results				Remarks ⁷				
											Xpert ⁴	Smear microscopy ⁶							
												1	2	3					

- 1 For diagnostic testing employing serial sputa or other specimens this is the date of receipt of the first set of specimens
- 2 Y=Yes; N=No; Unk = unknown
- 3 Y = previously treated; N = not previously treated, Unk = unknown
- 4 Patient on TB treatment; indicate months of treatment at which follow-up examination is performed
- 5 Xpert MTB/RIF test result reported as follows :
 - T = MTB detected, rif resistance not detected;
 - RR = MTB detected, rif resistance detected;
 - TI = MTB detected, rif resistance indeterminate
 - N = MTB not detected;
 - I = Invalid / no result / error
- 6 Smear results reported as follows:
 - 0=No AFB;
 - Scanty (and report number of AFB) = 1-9 AFB per 100HPF;
 - + = 10-99 AFB per 100 HPF;
 - ++ = 1-10 AFB per HPF;
 - +++ = > 10 AFB per HPF
- 7 If Xpert MTB/RIF indeterminate result, indicate error code or 'invalid'



The Laboratory Numbering system

The LN begins at number “1” at the start of each year.
It increases by one with each patient, until the end of the year.

Do not return to LN 1 at the end of each day, week, or month

Master List

Include a completed Master List whenever you send specimens to the smear microscopy laboratory.


Master List		
Laboratory name:		
Address:		
Patient name	Specimen 1	Specimen 2
Packed by		
<i>Name</i>	<i>Signature</i>	
<i>Dispatched Date</i> / /	<i>Time</i> : AM/PM	

Abbreviation	Name in full
AFB	Acid-Fast Bacilli
BSC	Biological Safety Cabinet
CDC	Centre for Disease Control and Prevention
CF	Carbol Fuchsin
EQA	External Quality Assessment
FM	Fluorescence Microscopy
GLI	Global Laboratory Initiative
JATA	Japan Anti-Tuberculosis Association
KNCV	KNCV Tuberculosis Foundation
KPIs	Key Performance Indicators
LED	Light Emitting Diode
LN	Laboratory Number
NTP	National Tuberculosis Programme
PPE	Personal Protective Equipment
QA	Quality Assurance
QC	Quality Control
RI	Refractive Index
TB	Tuberculosis
The Union	International Union Against Tuberculosis and Lung Disease
v/v	Volume for volume
VWS	Ventilated Work Station
WHO	World Health Organization
ZN	Ziehl-Neelsen

Laboratory design

The basic requirements for a sputum microscopy laboratory include:

1. Good ventilation
Directional ventilation provides healthy air for breathing. Air that may be contaminated by laboratory processes should flow away from staff and out of the laboratory.
2. A strong table/bench to prepare smears
3. A sink or plastic basin to stain smears
4. A table/bench to examine smears
5. A table/bench for paperwork
6. Basin for hand washing
7. Good lighting
8. Non-slip flooring
9. An area for receiving specimens



Biological Safety Cabinets

A biological safety cabinet (BSC) is not required for sputum smear microscopy.

- Only laboratories performing culture and drug susceptibility testing need a functioning BSC
- Never use a clean air cabinet, it can blow TB organisms into the laboratory

Contamination and infection control

Assume all samples are potentially infectious

Aerosols

Good work practice minimises aerosol formation and contamination of work surfaces and equipment. (See page 8 Personal Safety).

Disinfection and Spills

Disinfection

Disinfectants recommended for use in TB laboratories contain phenols, chlorine or alcohol.

Disinfection methods	Surfaces	Spills	Prepare
Phenol 5%	Yes	Yes	Every 2 days
Alcohol 70% v/v	Yes	No	Weekly
Hypochlorite 0.5%	No	Yes	Every 2 days

If your skin is contaminated with phenol, bleach or alcohol, wash thoroughly with soap and water

Phenol

- Toxic if swallowed
- Phenol is highly irritating to the skin, eyes and mucous membranes (e.g. lungs)
- Due to its toxicity and smell synthetic phenol derivatives are generally used in place of phenol

Chlorine

- Bleach is highly alkaline and will corrode metal
- Sodium hypochlorite solutions (domestic bleach) contain 35-150 g/l available chlorine – store in a well ventilated dark area
- Dilute in water to obtain a final concentration of 0.5%

Alcohol

- Volatile and flammable
- Keep away from open flames
- Store in proper containers to avoid evaporation
- Label bottles clearly – do not autoclave





Spills

Treat all spills as potentially infectious

1. Put on a laboratory coat and gloves
2. Place paper towel or cloth over the spill area and liberally apply disinfectant solution
3. Leave covered – minimum 15 minutes
4. Clean up the contaminated material and put into the waste container
5. Clean with a final wash using 70% v/v alcohol
6. Wash your hands after the clean up is complete



Waste management

Treat all laboratory waste as infectious

Laboratory staff are responsible for waste management and ensuring that anyone who must handle waste, including cleaners, drivers etc. is properly trained.

Where available autoclave laboratory waste before disposal.

Place potentially infectious waste into bins that have a disposable plastic lining with disinfectant added. When moving waste within or outside the laboratory, put it into a larger leak-proof plastic bag, tied at the top.


Laboratory staff are responsible for ensuring safe movement of laboratory waste.

When moving waste outside of the laboratory, the waste should be sealed in a container with a lockable lid.



Waste disposal

Locate the burning drum away from people in an open area as the fumes are toxic






Ergonomics

Good ergonomics reduces fatigue and injury



Good posture – supporting your feet straightens your back



Poor posture – feet unsupported



Good posture – raise the microscope to help straighten your back and keep your feet flat on the floor




Poor posture – seat too high or bench too low – feet not flat

Ventilated workstation

A ventilated workstation (VWS) is a partially enclosed workspace. Air is drawn inward, away from the technician and exhausted outside the laboratory, VWS are inexpensive to build and require little maintenance. VWS do not replace careful attention to risk minimising laboratory methods.

For more information on VWS see *Ventilated Workstation Manual for AFB Smear Microscopy* (see page 83).



Accurate laboratory results rely on internal monitoring (Quality Control and Key Performance Indicators) and EQA.

Why do Quality Control?

The purpose of Quality Control (QC) is to ensure that staining solutions work well and that they are not contaminated with AFB. Good quality solutions and staining technique make reading and reporting easier and more reliable. Accurate record keeping of preparation and testing provides confidence in your results.

Technicians preparing new staining solutions are responsible for QC before the solutions are used.

Technicians performing AFB-staining are responsible for regular QC using positive control smears.

Technicians who prepare control smears are responsible for their QC.

Preparing unstained control smears

Positive control smears

Ideal positive control smears are easy to count low-positives in the 1+ range.

1. Confirm a 1+ result for the selected specimen on 2 or 3 stained smears:
 - After liquefaction (standing overnight) and
 - Mixing – with sputum pot closed
2. Make at least 50 even equally sized smears from this confirmed 1+ sample, and air dry
3. Heat fix
4. On each slide write the positive control batch number and serial number within the batch
5. Check the number of AFB:
 - Randomly select six smears from this batch
 - Stain and carefully count the AFB
6. Start a separate page in your logbook for quality control of staining solutions
7. Record the batch number and results for each of the six smears then calculate the average number of AFB per smear length or per field
8. Store smears in a closed slide box labelled “*Positive control smears*”

Negative control smears

Make negative control smears from egg white diluted 5% in distilled water.

1. To assist focusing mix with a little sputum or saliva (containing cells)
2. After staining check a few smears to make sure there is no contamination with AFB
3. Make at least 50 even equally sized smears from this sample, and air dry
4. Heat fix
5. On each slide write the negative control batch number and serial number within the batch
6. Store smears in a closed slide box labelled “*Negative control smears*”

Testing solutions

To test the performance of a *freshly prepared solution* stain and examine:

- Two positive control smears *stained once* and
- Two negative control smears *stained three times*

The other solutions required can also be those already in routine use.

1. Stain negative smears *three times* to check for environmental mycobacteria
 - Only repeated staining makes these contaminants visible
2. Examine control smears carefully for:
 - Number and intensity of AFB colour
 - Complete de-colourisation of background
 - Absence of crystals and primary stain coloured artefacts
3. Compare your count with the number of AFB expected for the batch of positive control smears
 - There should be no negative or very low counts
 - AFB should show strong, solid colour
4. Accept the batch if it passes on all these points

Unsatisfactory results

1. Check the preparation technique, the quantities and reagents used:
 - If results are uncertain, stain a few more control smears making sure your technique is correct
2. Accept the batch if results are good
3. If results are again unsatisfactory:
 - Discard the bad batch of staining solution
 - Record the reason for rejection
 - Prepare fresh solution and perform QC

Keep accurate QC records in the logbook for all solutions prepared. Good records serve as an important reference to defend against possible complaints.

Monitoring

Key Performance Indicators (KPIs) are useful for internal and external evaluation of AFB-microscopy quality. They should be calculated monthly or quarterly from the Laboratory Register counts, and the results recorded in a chart.

Each laboratory is responsible for calculating its KPIs.

Monitoring trends within the laboratory should alert staff to identify a shift from normal patterns. Values that are too high or low may indicate a problem, however the acceptable range depends on the setting.

The TB Programme should collect individual laboratory KPIs and compare them across the laboratory network. This allows each laboratory to compare their performance with similar laboratories in the same area.

Reporting the data


	Negative	Scanty	1+	2+	3+	Total positive or scanty	Total smears
Suspect smears	a	b	c	d	e	$f=(b+c+d+e)$	$g=(a+f)$
Follow-up smears	h	i	j	k	l	$m=(i+j+k+l)$	$n=(h+m)$

Calculations

Workload	$g+n$
% positive suspect smears	f/g
% positive follow-up smears	m/n
% low positive suspect smears	$(b+c)/f$

Plot KPI's monthly or quarterly to obtain a trend line.

Plotting may not be effective if denominators (totals) are very small.



Target values

Laboratories should aim for:

- TB suspects – about 10% positives
- Follow-ups – about 5-10% positives
- Low positive suspect smears – about 30-50% of all positive suspect smears

EQA

EQA of AFB-microscopy commonly includes rechecking a randomly selected subset of routine smears by an external agency. For EQA to be effective technicians should keep all smears until the subset of smears has been selected and removed for rechecking.

The EQA process

When preparing slides for examination:

- Label all slides clearly with the LN and sample number
- Let oil soak into absorbent paper overnight after reading
- Store in numerical sequence leaving a space for the smear of the second sample
- Never write results on the slide

After the subset of routine smears has been selected for EQA and removed for rechecking, the remaining slides can be discarded.

Reuse the slide racks to start a new collection of routine slides. Store slides in numerical order leaving a space for the smear of the second sample.

Leave a space for second sample



Angra P, Becx-Bleumink M, Gilpin C, Joloba M, Jost K, Kam KM, Kim SJ, Lumb R, Mitarai S, Ramsay A, Ridderhof J, Rieder HL, Selvakumar N, van Beers S, van Cleeff M, Van Deun A, Vincent V. Ziehl-Neelsen staining: strong red on weak blue, or weak red under strong blue? *Int J Tuberc Lung Dis* 11: 1160-1161, 2007

Aziz MA, Ba F, Becx-Bleumink M, Bretzel G, Humes R, Iademarco MF, Kim SJ, Lamothe F, Paramasivan CN, Ridderhof J, Sloutsky A, Van Deun A, Shah KV, Weyer K. External quality assessment for AFB smear microscopy, Washington, DC:Association of Public Health Laboratories, 2002

CDC/The Union/GLI/APHL. Ventilated Workstation Manual for AFB Smear Microscopy: Manufacturing, Validation and User Guide. 2011
http://www.aphl.org/aphlprograms/global/Documents/GH_2011July_VentilatedWorkstationGuidance.pdf

International Union Against Tuberculosis and Lung Disease. Technical guide: Sputum Examination for Tuberculosis by Direct Microscopy in Low Income Countries. 2000

Van Deun A, Aung KJM, Hamid Salim A, Gumusboga M, Nandi P, Hossain MA. Methylene blue is a good background stain for tuberculosis light-emitting diode fluorescence microscopy. *Int J Tuberc Lung Dis* 14:1571-1575, 2010

Van Deun A, Hossain MA, Gumusboga M, Rieder HL. Ziehl-Neelsen staining: theory and practice. *Int J Tuberc Lung Dis* 12:108-110, 2008


World Health Organization. Laboratory services in tuberculosis control. Microscopy. Part II 1998. WHO/TB/98.258
http://www.who.int/entity/tb/publications/who_tb_98_258/en/index.html

World Health Organization. Policy statement: Same-day diagnosis of tuberculosis by microscopy. 2011. WHO/HTM/TB/2011.7
http://whqlibdoc.who.int/publications/2011/9789241501606_eng.pdf

World Health Organization. Policy statement: Fluorescent light-emitting diode (LED) microscopy for diagnosis of tuberculosis. 2011. WHO/HTM/TB/2011.8
http://whqlibdoc.who.int/publications/2011/9789241501613_eng.pdf

World Health Organization. Global tuberculosis report 2012. WHO/HTM/TB/2012.6
http://apps.who.int/iris/bitstream/10665/75938/1/9789241564502_eng.pdf

World Health Organization. Tuberculosis laboratory biosafety manual. 2012
WHO/HTM/TB/2012.11
http://apps.who.int/iris/bitstream/10665/77949/1/9789241504638_eng.pdf



- Your doctor/nurse has sent you to the laboratory because they suspect that you may have the symptoms of tuberculosis (TB)
- To diagnose TB two sputum specimens are needed and they will be collected:
 1. At first presentation
 2. Next morning before breakfast
- Collect specimens in the open air
- Good quality specimens from the lungs are required **not** saliva or nasal secretions
- Rinse your mouth out with bottled water if you have recently eaten, or if you have dentures (remove them first)

